
USER MANUAL


SIRIO

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From **380^{nm}**
wavelength To **980^{nm}**
wavelength

CONTACT INFO

✉ support@toplifeproject.com

🌐 www.toplifeproject.com

🌐 www.elioslamp.com

🌐 www.evolutamente.it

TABLE OF CONTENTS

ENGLISH

<i>Technical specifications</i>	1
<i>What's in the box</i>	3
<i>How it works</i>	4
<i>Precautions</i>	5
<i>Operating modes</i>	6
<i>Treatment Protocols</i>	8
<i>Contraindications</i>	18
<i>Maintenance and Care</i>	19
<i>Warranty and support</i>	20

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Before using the device, please carefully read the manual to know all the functions to ensure safe use.

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USER MANUAL

SIRIO

**MULTI-SPECTRUM LAMP
FOR AESTHETICS, BEAUTY,
AND WELL-BEING.**

1 TECHNICAL SPECIFICATIONS

Units: nm = nanometer (unit of light wavelength), W = watt (power), Hz = hertz (light pulses per second), mW/cm² = milliwatt per square centimeter (unit of light intensity on the skin)

TOTAL POWER	200W
NUMBER OF LEDS	200
STRUCTURE	Aluminium
DIMENSIONS	35 cm (height) × 21 cm (max width) × 17 cm (min width) × 2 cm (thickness)
WEIGHT	1,4 kg
POWER SUPPLY	External transformer (no battery)
TIMER	5–60 minutes (5-minute steps)
EMISSION MODE	Continuous + Pulsed (1, 10, 20, 40 Hz)
FLICKER	Flicker Free
EMF EMISSIONS	Low emissions
CERTIFICATIONS	CE
WARRANTY	2 years
ESTIMATED LED LIFE	50,000 hours

Emission spectrum

WAVELENGTH	COLOUR	N° LED	POWER	MODE
380 nm	Violet	6	6W	BLUE
405 nm	Blue-violet	6	6W	BLUE
460 nm	Blue / Cyan	6	6W	BLUE
590 nm	Yellow	36	36W	RED
630 nm	Red	37	37W	RED
670 nm	Deep red	37	37W	RED
810 nm	Near-infrared (NIR)	36	36W	RED
980 nm	Near-infrared (NIR)	36	36W	RED


Emission spectrum — Irradiance

DISTANCE	RED + NIR	BLUE	RED+BLUE
10 cm	100 mW/cm ²	—	—
15 cm	93 mW/cm ²	—	—
20 cm	86 mW/cm ²	—	—
30 cm	70 mW/cm ²	31 mW/cm ²	90 mW/cm ²


Note: For BLUE and RED+BLUE mode the recommended protocol is at 30 cm — the optimal distance for the photobiological safety of the blue and violet wavelengths.

2 WHAT'S IN THE BOX


1. MAIN DEVICE

2. POWER ADAPTER

3. PROTECTIVE GOGGLES
(WRAPAROUND, ANTI-UV)

3 HOW IT WORKS

EN

SIRIO is a multifunction LED aesthetic lamp **designed for the wellness of the skin, eyes, hair and body.**

Thanks to eight precisely selected wavelengths — from violet to near-infrared — SIRIO **offers three distinct usage modes**, each designed for a specific purpose.

It is not a medical device. It is a personal wellness tool, designed to complement your daily care and regeneration routine.

The three emission modes

- **RED mode** — red, deep red and near-infrared light (+ yellow)
- **BLUE mode** — blue, blue-violet and violet light
- **RED+BLUE mode** — full spectrum, all wavelengths active

4 PRECAUTIONS

EYE PROTECTION

Always wear the protective goggles included in the package during every session. Never look directly at the LEDs, especially in BLUE mode (violet and blue are the most energetic wavelengths).

Usage timing — IMPORTANT

- **RED mode:** can be used at any time of day, including the evening.
- **BLUE and RED+BLUE modes:** recommended EXCLUSIVELY in the morning, before 13:00. Blue, cyan and violet wavelengths (380–460 nm) activate non-visual opsins and inhibit melatonin production. Using them from the afternoon onwards can disturb sleep and alter the circadian rhythm.


Simple rule: BLUE and RED+BLUE in the morning only. RED at any time.

Minimum safety distances

- RED + NIR mode: minimum 10 cm from the skin
- BLUE mode: minimum 30 cm from the skin
- RED+BLUE mode: minimum 30 cm from the skin

5 OPERATING MODES

EN

The SIRIO lamp has 4 buttons on the control panel:

- 1. RED:** activates red, deep red and near-infrared light mode (+ yellow)
- 2. BLUE:** activates blue, blue-violet and violet light mode
- 3. PULSE:** selects the pulse frequency (each press: OFF > 1 Hz > 10 Hz > 20 Hz > 40 Hz > OFF)
- 4. TIME:** sets the timer (each press adds 5 minutes, from 5 to 60 minutes)

At the centre of the panel is the timer display, which shows the selected session time and the remaining time.

Instructions for use — step by step

Before the session

- Place SIRIO on the table and tilt the panel to the desired angle
- Connect the transformer to the power outlet
- Put on the protective goggles when the blue lights are on

- Make sure the area to be treated is clean and free of make-up or creams

During the session

- Press RED, BLUE or both to select the desired mode
- Press PULSE to select the pulse frequency (or leave on Continuous)
- Press TIME to set the timer (each press = +5 minutes)
- Position yourself at the correct distance for the chosen mode
- Relax during the session — there is no need to move

After the session

- The lamp switches off automatically when the timer expires
- Remove the protective goggles
- You can apply your preferred moisturiser or serum (the mild warmth aids absorption)

6 TREATMENT PROTOCOLS

SECTION A — RED MODE (RED + NIR)

A1. What RED mode is

RED mode activates the wavelengths 590 nm (yellow), 630 nm, 670 nm (red), 810 nm and 980 nm (near-infrared — NIR). These wavelengths penetrate deep into the skin and underlying tissues, stimulating the mitochondria — the energy powerhouses of the cells — to increase energy (ATP) production. The result is faster cellular regeneration, reduced inflammation and improved microcirculation.

Yellow (590 nm) plays a special role: it acts on the surface capillaries and cutaneous microcirculation, and is particularly indicated for rosacea and eye health.

A2. What RED mode is for

- Skin rejuvenation and aesthetic care — stimulates the production of collagen and elastin
- Microcirculation — yellow (590 nm) improves circulation in the surface capillaries
- Rosacea — yellow (590 nm) acts on dilated capillaries, reducing redness and couperose
- Eye health — indirect application (never direct): supports retinal health
- Hair health — stimulates scalp microcirculation and hair follicles
- Wounds and minor skin inflammation — accelerates healing

A3. Pulse frequencies — RED mode

FREQUENCY	BIOLOGICAL EFFECT	WHEN TO USE IT
Continuous (OFF)	Constant action on surface and deep tissues	General aesthetic treatments, rejuvenation, microcirculation, rosacea
1 Hz	Slow rhythm — promotes deep relaxation and cellular regeneration	Wounds, healing, chronic inflammation, evening treatments
10 Hz	Medium rhythm — balance between stimulation and regeneration	Hair loss, hair follicles, deep anti-aging
20 Hz	Sustained rhythm — cellular metabolic stimulation	Skin firmness, collagen, energising morning treatments
40 Hz	Fast rhythm — neuromuscular stimulation (gamma)	Active microcirculation, swelling reduction, lymphatic drainage

SECTION B — BLUE MODE (BLUE + VIOLET + CYAN)

Recommended timing: morning only, before 13:00. Blue, cyan and violet wavelengths activate non-visual opsins and inhibit melatonin. Using this mode from the afternoon onwards can disturb sleep and alter the circadian rhythm.

B1. What BLUE mode is

BLUE mode activates the wavelengths 380 nm (violet), 405 nm (blue-violet) and 460 nm (blue/cyan). Violet (380 nm) is at the boundary between visible light and UV-A. It has natural antimicrobial properties and acts on melanin and the skin's surface chromophores.

Cyan (460 nm) is the most powerful wavelength for regulating the circadian rhythm: it activates melanopsin in the retinal ganglion cells and in the cutaneous opsins, directly influencing melatonin levels and the sleep-wake cycle.

Protocol always at 30 cm from the skin. Irradiance at 30 cm: 31 mW/cm².

B2. What BLUE mode is for

- Deep facial detox — natural antimicrobial properties that purify the skin in depth
- Circadian reset — cyan (460 nm) and violet (380 nm) activate non-visual opsins, regulating the internal biological clock
- Acne — blue-violet (405 nm) inactivates *Cutibacterium acnes* bacteria through the production of photoactive porphyrins

B3. Pulse frequencies — BLUE mode

FREQUENCY	BIOLOGICAL EFFECT	WHEN TO USE IT
Continuous (OFF)	Constant action — maximum exposure to blue wavelengths	Acne treatment, deep detox, skin purification
1 Hz	Slow rhythm — gentle, prolonged action	Sensitive skin, light purification
10 Hz	Medium rhythm — balancing and regulation	Morning circadian reset, biological balance
20 Hz	Sustained rhythm — stimulation of non-visual opsins	Active circadian synchronisation, alertness
40 Hz	Fast rhythm — gamma stimulation	Cognitive activation, morning waking phase

SECTION C — RED+BLUE MODE (FULL SPECTRUM)

Recommended timing: morning only, before 13:00. The presence of blue, cyan and violet wavelengths (380–460 nm) inhibits melatonin production. For afternoon and evening treatments, use RED mode exclusively.

C1. What RED+BLUE mode is

RED+BLUE mode activates all 8 wavelengths simultaneously for a total of 200 W: from violet (380 nm) to deep red and near-infrared (980 nm). It is the most complete treatment SIRIO can offer. It combines the regenerating and anti-

inflammatory action of red and NIR with the purifying, circadian and antimicrobial action of blue and violet.

Protocol always at 30 cm. Total irradiance at 30 cm: 90 mW/cm².

C2. What RED+BLUE mode is for

- Complete aesthetic facial treatment — regeneration + purification + microcirculation in a single session
- Complete circadian reset — blue/cyan wavelengths regulate the biological clock while red and NIR regenerate tissues
- Hair health — advanced treatment: red/NIR stimulates the follicles, blue/violet purifies the scalp
- Wounds and minor inflammation — advanced treatment: synergistic regenerative + antimicrobial action
- General wellness — a combined session to make the most of all the benefits in a single sitting

C3. Pulse frequencies — RED+BLUE mode

FREQUENCY	BIOLOGICAL EFFECT	WHEN TO USE IT
Continuous (OFF)	Maximum action on all tissues	Complete aesthetic treatment, general wellness
1 Hz	Deep regeneration and relaxation	Skin recovery, wounds and inflammation
10 Hz	Biological balance — follicular stimulation + purification	Advanced hair care
20 Hz	Active stimulation — collagen + non-visual opsins	Complete morning anti-aging treatment

FREQUENCY	BIOLOGICAL EFFECT	WHEN TO USE IT
40 Hz	Maximum metabolic and circadian stimulation	Energising sessions, physical and mental preparation

SUMMARY TABLE THE 4 PULSE FREQUENCIES

FREQUENCY	WHAT IT DOES IN THE BODY	MAIN APPLICATIONS
Continuous	Constant, stable action. Maximum energy delivered to tissues.	General aesthetic treatments, acne, detox, microcirculation, rosacea
1 Hz — Delta waves	Ultra-slow rhythm. Promotes deep cellular regeneration. Similar to brain rhythm during deep sleep.	Wounds, healing, inflammation, evening treatments, delicate skin
10 Hz — Alpha waves	Medium rhythm. Balance between stimulation and regeneration. Typical of the brain in relaxed wakefulness.	Hair loss, hair follicles, deep anti-aging, balancing treatments
20 Hz — Beta waves	Sustained rhythm. Active metabolic stimulation. Optimal for the synthesis of collagen and elastin.	Anti-aging, skin firmness, collagen, morning treatments
40 Hz — Gamma waves	Fast rhythm. Maximum neurological efficiency and microcirculation stimulation.	Microcirculation, swelling reduction, drainage, active circadian reset

Pulse frequencies indicate how many times per second the light switches on and off. Each frequency interacts differently with the body's biological systems.

How to choose: always start with Continuous mode for general treatments. Use the pulse frequencies to enhance a specific effect: deep regeneration with 1 Hz, active stimulation with 20–40 Hz.

USAGE PROTOCOLS — RED MODE

Face and neck — Rejuvenation and anti-aging

Distance: 15–20 cm

Time: 10–15 minutes

Frequency: Continuous or 20 Hz

Sessions: 4–5 times a week for the first 30 days, then 2–3 times a week

Rosacea and couperose — Yellow (590 nm)

Distance: 20 cm

Time: 8–10 minutes

Frequency: Continuous or 1 Hz

Sessions: 3 times a week

Note: RED mode includes yellow (590 nm), which acts directly on the dilated capillaries typical of rosacea

Eye health

Distance: 30 cm

Time: 5 minutes

Frequency: Continuous or 10 Hz

Sessions: 3–4 times a week

IMPORTANT: The treatment is INDIRECT. Position the lamp to the side or from below, with eyes closed and protective goggles on.

Hair health — Scalp

Distance: 10–15 cm

Time: 10–15 minutes

Frequency: 10 Hz or 20 Hz

Sessions: 3–4 times a week

Note: Position the lamp above the scalp, on dry, clean hair.

Wounds and minor skin inflammation

Distance: 15–20 cm

Time: 5–8 minutes

Frequency: 1 Hz or continuous

Sessions: once a day until improvement

Only for small, superficial, non-infected wounds. In case of deep wounds or wounds that do not improve, consult a doctor.

USAGE PROTOCOLS — BLUE MODE

Deep facial detox

Distance: 30 cm

Time: 10–15 minutes

Frequency: Continuous or 1 Hz

Sessions: 2–3 times a week

Note: Clean face, without make-up or creams. After the session, apply your moisturiser.

Morning circadian reset

Distance: 30 cm

Time: 10–20 minutes

Frequency: 20 Hz or 40 Hz

Sessions: every morning on waking

Note: Cyan (460 nm) synchronises the biological clock,

improving alertness, mood and the sleep-wake cycle.

Acne treatment

Distance: 30 cm

Time: 10–15 minutes

Frequency: Continuous

Sessions: 4–5 times a week for 4–8 weeks

Note: Clean face. Do not apply photosensitising cosmetic products before the session.

USAGE PROTOCOLS RED+BLUE MODE

Complete aesthetic facial treatment

Distance: 30 cm

Time: 15–20 minutes

Frequency: Continuous or 20 Hz

Sessions: 3–4 times a week

Hair health — Advanced treatment

Distance: 30 cm

Time: 15 minutes

Frequency: 10 Hz

Sessions: 3–4 times a week

Complete circadian reset (morning)

Distance: 30 cm

Time: 15–20 minutes

Frequency: 20 Hz or 40 Hz

Sessions: every morning

Wounds and minor inflammation — Advanced treatment

Distance: 30 cm

Time: 5–8 minutes

Frequency: 1 Hz

Sessions: once a day (morning only)

Note: Red/NIR accelerates cellular regeneration, blue/violet exerts an antimicrobial action. Only for small, superficial, non-infected wounds. In case of deep wounds or wounds that do not improve, consult a doctor.

7 CONTRAINDICATIONS

EN

Who should NOT use SIRIO without medical advice:

- Pregnant women
- Children under 12
- People with photosensitive epilepsy
- People with confirmed or suspected skin tumours
- People taking photosensitising medication (antibiotics, antidepressants, retinoids — consult a doctor)
- People with implantable electronic devices (pacemakers, defibrillators)
- People with serious eye conditions or retinopathies

Do NOT use SIRIO:

- On deep open skin lesions or recent surgical wounds
- Directly on the eyes without protective goggles
- On areas with recent tattoos (less than 4 weeks)
- In combination with photosensitising products applied to the skin

8 MAINTENANCE AND CARE

- Clean the surface of the LED panel with a soft, slightly damp cloth
- Do not use aggressive detergents, alcohol or spray directly on the panel
- Never immerse the lamp in water
- Store the protective goggles in a safe place
- Never disassemble the device — in case of malfunction, contact Elioslamp support

STORAGE AND DISPOSAL

- Store SIRIO in a dry place, away from heat sources and humidity
- Do not expose to extreme temperatures (below 0 °C or above 40 °C)
- Disposal: SIRIO is an electrical appliance — do not dispose of it with ordinary waste. Take it to the authorised WEEE collection centres in your municipality.

9 WARRANTY AND SUPPORT

EN

Dear Customers,
Thank you for choosing our product. Customer satisfaction is our primary commitment. For your peace of mind, we offer:

- A limited 2-year warranty from the date of purchase
- Warranty coverage begins from the purchase date indicated on your invoice
- Please keep the original invoice for any warranty claims

WARRANTY EXCLUSIONS

The warranty does not cover:

- Damage caused by improper operation or use.
Malfunctions caused by disassembly of the device.
- Malfunctions or damage caused by accidents or human activities (scratches, immersion in water, or unauthorized replacement of parts).
- Damage caused by natural disasters.

For technical issues, malfunctions, or warranty claims, please contact:

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